

# PsyMSc4 Kog - Praxismodul

# Python for Psychologists

```
def dotwrite(ast):  
    nodename = getNodeName(ast[0],ast[0])  
    label=symbol.symbol(nodename,ast[0])  
    print '%s [%s] (%s)' % (label,ast[0],label),  
    if isinstance(ast[1],list):  
        if ast[1]:  
            print '  
        else:  
            print ''  
    else:  
        print ''  
        children = []  
        for n, child in enumerate(ast[1:]):  
            children.append(dotwrite(child))  
        print '%s-> [%s]' % (nodename,children),  
        for name in children:  
            print '%s' % name,
```



python

Rebecca Mayer  
Dominik Kraft  
WS 19/20




**... a module for programming experiments**

# PsychoPy - Getting started I

Builder View

Coder View



# PsychoPy - Where to get help

some useful resources:

- **PsychoPy Manual** <https://www.psychopy.org/PsychoPyManual.pdf>
- **API Reference Manual** <https://www.psychopy.org/api/api.html>
  - Contents:
    - `psychopy.core` - basic functions (clocks etc.)
    - `psychopy.visual` - many visual stimuli
    - `psychopy.clock` - Clocks and timers
    - `psychopy.data` - functions for storing/saving/analysing data

# Experiment Header

```
>> #!/usr/bin/env python
```

→ tells your OS that this programme is using the python language

```
>> # -*- coding: utf-8 -*-
```

→ character encoding

```
# Import the PsychoPy libraries that you want to use
```

```
>> from psychopy import core, visual, event
```

## Presenting stimuli

`psychopy.visual.window(size=[800,600], pos=None, color=(0,0,0), colorSpace='rgb', ...)`


- creates a window in which stimuli etc. can be presented
- important method: `.flip()`
  - `.flip()` returns the time of presentation

## Presenting stimuli

`psychopy.visual.TextStim(window, text="Hello world", color=(0,0,0), colorSpace='rgb', pos=(0,0), ...)`


- creates a text object
- important method: `.draw()`

# Presenting stimuli







# Presenting stimuli



# Presenting stimuli



# Presenting stimuli



# Example

```
>> from psychopy import core, visual
```

```
# Create a window
```

```
>> win = visual.Window([400,300], monitor="testMonitor")
```

```
# Create a text stimulus for a certain window
```

```
>> message = visual.TextStim(win, text="Hi there!")
```

```
# Draw the stimulus to the window.
```

```
>> message.draw()
```

```
# Flip backside of the window. Puts picture immediately on the screen (but only for a fraction of time)
```

```
>> win.flip()
```

```
# Pause 5 s, so you get a chance to see it!
```

```
>> core.wait(5.0)
```

```
# Finally, close the window
```

```
>> win.close()
```

## Timing stimuli

# tells psychopy to wait t seconds, before doing anything else

```
>> core.wait(t)
```

# initialize a clock for your experiment

```
>> clocki = core.clock()
```

# returning the current time of your clock

```
>> clocki.getTime()
```

# reset your clock (we named it clocki)

```
>> clocki.reset()
```

## Recording responses

`psychopy.event.waitKeys(maxWait=inf, keyList=None, timeStamped=False, ...)`

- waits maxWait seconds for either any key to be pressed or for one of the keys contained in keyList to be pressed (if specified)
  - (Moves on, if the key is pressed before maxWait is reached)
- timeStamped
  - False: only the key pressed is returned (as a list); if no key is pressed, None is returned
  - True: the key pressed is returned as well as the time (as a tuple) at which the key press occurred; again, None is returned, if no key was pressed
  - Clock: give the name of your clock to have the clock time saved instead; again, None is returned if no key was pressed

## Exercise

- 1) Import the libraries core, visual and event from the psychopy module
- 2) Create a window of 500 x 500. Also, close the window again at the end of your script.
- 3) Now, create a text stimulus (a letter of your choice) that is presented on the screen (use the default settings).
- 4) Use `.waitKeys()` and the `maxWait` argument to a) record the response and, thus b) show the stimulus for 3 seconds. For now, don't specify any `keyList` and set `timeStamped` to `False`. Also, save the key to a variable (fyi: the output of `.waitKeys()` is of length 1). Note: you should also cover the case of misses, otherwise you might get an error when trying to assign the output. Hint: If no key has been pressed you the return value is `None`.
- 5) Now, as a second step, also record the reaction time. As a first step, simply set `timeStamped` to `True`. Store both the key and the reaction time in separate variables. Note, that now the output of `wait.Keys()` is of length 2, so you will have to "offer" a tuple of two variables, which the output can be assigned to.
- 6) Finally, write a loop. In each iteration, have another letter presented to the screen (check the `setText` method). Also, store all keys and rts in a list.
- 7) Save the presentation times of your stimuli (that is, at what time they have been presented).

## Timing stimuli

# tells psychopy to wait t seconds, before doing anything else

```
>> core.wait(t)
```

# initialize a clock for your experiment, your timer starts right away!

```
>> clocki = core.clock()
```

# returning the current time of your clock

```
>> clocki.getTime()
```

# reset your clock (we named it clocki)


```
>> clocki.reset()
```



# The Spyder IDE

variable explorer/...

editor



console